carilex[®]

VT · 200-i

Negative Pressure Wound Therapy System
Medical Powered Suction Pump with Canister
Instructions for Use

About This Document

Congratulations and thank you for purchasing this high quality NPWT. Please read these Instructions carefully before use and observe the safety instructions and requirements for the proper operation and maintenance of the device.

Identification Details of the Device

Device identification

These Instructions for Use are intended exclusively for devices with the following specification:

Device name: Negative Pressure Wound Therapy System - Medical Powered Suction Pump with Canister

The serial number is shown on the label on the rear of the power unit.

Details of the Device Documentation

Validity of the documentation

This manual describes the VT · 200-i Powered Suction Pump. It is part of the device documentation. Do not share this device with a third party without including these Instructions for Use.

The end page of these Instructions for Use shows the current document edition.

Subject to change

The contents of the Instructions for Use can be changed by the manufacturer at any time without prior notice.

Translations

The English version of these Instructions for Use is authentic. In the event of any clarification, question or dispute as to the content of any translation of these Instructions for Use, the English version shall prevail.

Copyright

The resulting rights for reprinting, photomechanical, digital processing, or copying of the content, or any part thereof, are permitted only with prior written approval from Carilex Medical Inc.

This provision shall not affect the reproduction for internal use.

◆ About This Document

Manufacturer

Carilex Medical, Inc.
No. 77, Keji 1st Rd., Guishan Dist., Taoyuan City 333, Taiwan

EU Representative:

Emergo Europe B.V.
Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands

If the user or/and patient occur any serious incident in relation to the device should be reported to Carilex and the local Authorized Representative immediately.

Sales and
service

For Support or Complaints

If you have any questions or concerns about the device or need accessories, please contact the distributor that delivered the device to you or your patients.

Table of Contents

1	Important Notes for Safe Use	5-14
1.1	Designation of the Groups of Individuals	5
1.2	Notes for the Users	6
1.2.1	Instruction	6
1.2.2	Transferring the Device	6
1.2.3	Maintenance and Installation	6
1.2.4	Information and Test Obligation of the User	6
1.3	Procedures for Accident Prevention	6
1.3.1	Procedures for Transferring the Device	6
1.3.2	Qualification Requirements Hygiene Staff	7
1.3.3	Availability of the Instructions for Use	7
1.4	Purpose of the Device	7
1.4.1	Intended Purpose	7
1.4.2	Precautions	8
1.4.3	General Precautions for Indication for Use	9
1.4.4	Caution	9
1.4.5	Exclusion Clause for Use	10
1.4.6	EMC Caution Statement	10
1.4.7	Safety Caution Statement	11
1.5	Warranty	12
1.6	Safety Information and Symbols of These Instructions for Use	12
1.6.1	Identification of Risks of Injury	12
1.6.2	Identification of Material Damage	12
1.6.3	Identification of Additional Information	12
1.6.4	Additional Symbols to the Safety Information	13
1.7	Safety Instructions before Use	13
2	Delivery and Storage	15
2.1	Packaging	15
2.2	Models	15
2.3	Delivery Control	15
2.4	Operating & Transport & Storage Conditions	15
3	Device and Functional Description	16-18
3.1	Device Description	16
3.2	Functional Description	17
3.2.1	Panel of the Power Unit	17
3.2.2	Function of the Power Unit	17
3.2.3	Therapy Mode	17
3.2.4	Symbols on the Display	18
4	Getting Started	18-21
4.1	Battery	18
4.2	Collection Canister	19
4.3	Dressing	20
4.4	Lay Operator Briefing Information	21

◆ Table of Contents

5	Operation Procedure	22-28
5.1	Check Points before Using Powered Suction Pump	22
5.2	Insert Canister into the Powered Suction Pump	22
5.3	Connect the Powered Suction Pump Tubing Connector to the Carilex VT Dressing Kit Tubing Connector, and Screw Tightly to Ensure a Tight Airlock.	22
5.4	Using the Function Keys	23
5.5	Turn the Unit On / Off	23
5.6	Therapy Mode Setting	23
5.7	Starting Negative Pressure Wound Therapy (NPWT)	25
5.8	Mute / Return	28
5.9	Change Powered Suction Pump Canister	28
6	Error Indicators	29
6.1	Error Indicators	29
7	Application of the Negative Pressure Wound Therapy	30
7.1	Application of the Negative Pressure Wound Therapy	30
8	Disinfection and Cleaning	31-32
8.1	Disinfection and Cleaning	32
9	Care and Maintenance	33
9.1	Inspection	33
9.2	Maintenance	33
10	Troubleshooting	34
11	Technical Data	35
12	EMC Declaration	36-37
13	Recommendations for Li-Ion Battery Pack	38

1 Important Notes for Safe Use

Named groups
of persons

1.1 Designation of the Groups of Individuals

The named groups of persons in these Instructions for Use are as follows.

Operators

An operator (surgical supplier, health insurance, clinic, etc.) is any legal person who owns and uses a Powered Suction Pump, or on whose behalf the device is used. The operator is responsible for providing a safe device and to instruct users properly on the operation and safe use of the device.

Users

Users are people who have been trained to:

- Operate the Powered Suction Pump.
- Supervise patients using the device for therapy or care purposes.

Users are fully responsible for the safe and correct use of the device. A review of the functions must be carried out and the proper conditions of the device must be confirmed by the user before each use or transfer for use.

Professionals

These authorized persons are skilled personnel who may be employees of the operator and who:

- Have acquired their knowledge through professional training in the medical-technical field,
- Carry out their activity on the basis of professional work experience and instructions according to safety-related regulations.
- Are able to detect possible hazards during work.

In countries in which the pursuit of an activity in the medical-technical area is certified, classification of qualified personnel is subject to appropriate approvals.

Patients

Patients are persons in need of care and who use the Powered Suction Pump for therapy or care purposes.

Lay Operator

A lay operator is a person without relevant specialized training. The lay operator should contact the manufacturer or manufacturer's representative under following conditions:

- For assistance in setting up, using or maintaining the ME Equipment System; or
- To report unexpected operation or events.

1 Important Notes for Safe Use

1.2 Notes for the Users

Note that the medical device should only be used by persons who have been trained in the operation and the intended purpose of the device.

Training on the device

1.2.1 Instruction

Training on how to use the device must be carried out by qualified personnel. If the device is approved for use, then users must comply with the Instructions for Use.

Device approval

1.2.2 Transferring the Device

The device may be used only if the operator has released it for therapeutic or care use and if the transfer was carried out properly and under the supervision of authorized personnel.

After the transfer, users are fully responsible for the safe and dedicated use of the device.

1.2.3 Maintenance and Installation

The maintenance and/or repair of the equipment or parts must be carried out only by an authorized service agent.

Obligation of the user

1.2.4 Information and Test Obligation of the User®

Read these Instructions for Use carefully before the first use of the device. This will allow you to experience all the benefits that the device offers and avoid possible personal injury and property damage.

A review of the functions must be carried out and the proper condition of the device must be confirmed by the user before each use or transfer for use by patients.

In case of specific issues that are not covered in enough detail in these Instructions for Use, please contact the distributor or operator for further guidance.

1.3 Procedures for Accident Prevention

The Powered Suction Pump is made according to current, reliable, state of the art technology. However, hazards may arise if it is operated by untrained personnel or it is not operated as described in these Instructions for Use.

1.3.1 Procedures for Transferring the Device

In order to comply with the regulations of accident prevention and to prevent accidental damage, follow these procedural guidelines when transferring the device:

- The device must be thoroughly cleaned and disinfected before the first use.
- The initial start-up of the device, as well as the transfer to the user, must be carried out by authorized personnel assigned by the operator.
- After completion of the training, it must be documented that the user understands the operation and use of the device for therapy or care purposes.

1 Important Notes for Safe Use

Hygiene measures

1.3.2 Qualification Requirements Hygiene Staff

The nature of hygiene measures is determined by the use environment of the device.

- If the device is used in clinical areas (e.g. hospitals, clinics, nursing homes, elderly homes, etc.) cleaning and disinfection must be carried out on the product or parts only by appropriately qualified personnel who are familiar with relevant hygiene regulations.
- When using the device in non-clinical areas, users or trained cleaning personnel may clean the device.

Obligation to provide information

1.3.3 Availability of the Instructions for Use

The Instructions for Use are an integral part of the device and must be stored in a place so that the safety instructions and other important information are accessible at any time and can be reviewed by users.

Do not pass the device to a third party without these Instructions for Use. Using the edition version as a guide always ensures that the current and valid Instructions for Use document is supplied with the device.

1.4 Purpose of the Device

To ensure the safety of patients and users, the device must be used only for its intended purpose.

1.4.1 Intended Purpose

Carilex "Negative Pressure Wound Therapy System" (NPWT) is indicated for patients who would benefit from wound management via the application of negative pressure for removal of fluids and excess exudates, irrigation fluids, infectious material, and tissue debris which may promote wound healing. To help provide safe and effective use, Carilex pump is to be used only with the Carilex disposable components: canister and dressing kits.

Intended users: Medical practitioner (Clinician use only)

Intended patient populations: No restriction in patient population.

Use environment: Carilex "NPWT" is intended to be used in hospitals, clinics, nursing home, elderly homes.

Indications:

- Traumatic
- Dehisced wounds
- Partial thickness burns
- Chronic wounds such as pressure ulcers, diabetic foot ulcers, venous leg ulcers
- Acute wounds
- Flaps and grafts.

Contraindications:

Patients with the following conditions:

- Presence of necrotic tissue
- Malignancy (except for quality of life reason for terminal patients)
- Exposed arteries, veins, nerves, organs or vascular graft
- Use over anastomotic sites

1 Important Notes for Safe Use

- Untreated osteomyelitis
- Unexplored or non-enteric fistulas

Relative Contraindications

- Ischemic wounds
- Ongoing infection
- Fragile skin
- Adhesive allergy

Clinical benefit:

NPWT benefits include rapid wound healing, reduction of dressing changes, reduced infection risk, reduced treatment costs, control of exudate, reduction of oedema and provision of a closed moist wound healing environment, concurrent rehabilitation, and better patient comfort and tolerance.

1.4.2 Precautions

The following statements describe medical conditions that may require special care to be exercised by a practitioner for the safe and effective use of the Powered Suction Pump.

1. Difficult wound hemostasis, or who are on anticoagulants
Patients on anticoagulation medicine or who have active bleeding or who have difficult wound hemostasis should be treated with caution. These patients are at an increasing risk for bleeding and bleeding complications and should be treated and monitored by properly trained medical caregivers in a controlled setting.
2. Exposed tendon, nerves or blood vessels should be protected
Close proximity of blood vessels, organs, muscle, and fascia. All blood vessels, organs, muscles, and fascia that are in close proximity to the wound site and/or are exposed and/or are near the skin surface should be properly protected prior to initiating NPWT. Patient with infections in the wound and or other parts of the body have to receive proper systemic treatment.
3. Weakened, irradiated or sutured blood vessels or organs
These patients are at an increasing risk for bleeding and bleeding complications and should be treated and monitored by properly trained medical caregivers in a controlled setting.
4. Bone fragments or sharp edges
Sharp edges from bony fragment may puncture blood vessels, organs, muscles, and fascia and may lead to bleeding. Proper care should be taken to cover the bony fragments and protect the wound area and other areas from bleeding.
5. Infected wound
Patient with infections in the wound and or other parts of the body have to receive proper systemic treatment. Infected wounds may need more frequent dressing changes, up to twice a day, and the patient and wound must be inspected regularly for signs of increased infection or sepsis.
6. For patients a known history of autonomic dysreflexia, please increase number of monitoring during the treatment as well as inspection for displacement of dressings.
7. Do not use NPWT if person experiences autonomic dysreflexia.

1 Important Notes for Safe Use

1.4.3 General Precautions for Indication for Use

1. It is important that a physician or other qualified healthcare provider evaluates the patient to ensure that the use of the Powered Suction Pump is an appropriate therapy.
2. To reduce the risk of transmission of blood-borne pathogens, regardless of their diagnosis or presumed infection status, all users should take medical standard operating procedure precautions against infection control.
3. Caregivers should wear gloves, a gown, and goggles if there is the possibility of contact with the patient's body fluids.
4. Change the dressing if the pump has stopped for more than two hours.
5. Consider mode of therapy (continuous, intermittent and instill mode).

NOTE

This product is for use only by individuals who have been adequately trained in using NPWT devices and who have had medical training in wound care. Operating this device or changing the settings should be done only by a physician's order or other qualified clinical caregiver.

1.4.4 Caution

The following caution statements describe the potential for serious consequences to the patient such as death, injury, or adverse reactions. Failure to read and follow all instructions in this manual prior to use may result in death or injury of the patient.

1. Physician should consider the patients' size and weight when prescribing this device. Infants, children, certain small adults and elderly patients should be closely monitored for fluid loss and dehydration.
2. The device is not safe for use with an MRI or PET scan and must be disconnected from the patient prior to MRI or PET scan.
3. The device is not safe for use with a Hyperbaric Oxygen Therapy (HBO) chamber and must be disconnected from the patient prior to entering the chamber.
4. The device may be used in the event that defibrillation is needed, provided there is no electrical connection between the patient and the device. In such case, the device must be completely disconnected from the patient. Be especially vigilant about removing wound dressing if defibrillation is required in the area of dressing placement. Failure to remove the dressing may inhibit transmission of electrical energy and/or patient resuscitation.
5. Residues of gauze/foam that may increase the risk associated with wound infection and bleeding. To prevent unintentional gauze/foam retention, all dressings should be carefully removed from the wound and the entire wound bed. Upon removal of the dressings, the wound bed should be cleaned in accordance with standard wound care practices (or facility guidelines), prior to the application of new sterile dressing.

1 Important Notes for Safe Use

6. If necessary, all wounds should be debrided prior to application of the therapy and/or dressings.
7. Ensure that there are no pockets left in the wound after application of the dressings.
8. Apply the dressing on the exposed arteries, veins, nerves, organs, or anastomotic site, that may increase the risk associated with wound bleeding. Please properly protect the wound site prior to initiating NPWT.
9. Apply the device over the unexplored an ischemia of soft tissue that may increase the risk associated with ischemic fasciitis. Please check the contraindications and properly protect the wound site prior to initiating NPWT.
10. Dressing does not be sealed that may increase the risk associated with wound infection. Please put the pump carefully in the carrying bag or fixed it, to avoid that excessively tension act on tube after pump fall and then leading the dressing come off from the skin.
11. Strangulation resulting from canister or dressing tube.
12. Skin irritation due to prolonged exposure to APPLIED PARTS or other ACCESSORIES.
13. To minimize the risk bradycardia, this device should not be placed in proximity to the vagus nerve.
14. Care should be taken if a patient has a spinal cord injury (potential for stimulation of the sympathetic nervous system).
15. Use of NPWT presents a risk of tissue ingrowth; users should be informed of means to reduce this risk (i.e. reducing therapy pressure, increasing the frequency of dressing changes, and monitoring).
16. This Powered Suction Pump is Not AP / APG protected. Do not use it under Oxygen enriched environment and flammable anesthetics.
17. Unapproved stop the therapy risks to patients because they have not been evaluated by the physician for wound assessment. Without physician evaluate, there is no way to know if the wound is healthy, whether the wound is infection or worsened. Do not stop the therapy by yourself, or, contact your caregiver if any problem.

1.4.5 Exclusion Clause for Use

Any and all applications outside of the conditions specified above is regarded as improper. The user and the operator respectively are exclusively liable for any damage resulting from improper use of the device.

1.4.6 EMC Caution Statement

1. Warning: Use of the VT · 200-i adjacent to or stacked with other RF (radio frequency) communication equipment (including antennas) should be avoided and used no closer than 30cm to any part of the VT · 200-i, including cables specified by the manufacturer. This could result in improper operation and/or degradation of the essential performance of this equipment.

1 Important Notes for Safe Use

2. The compliance statement of emission class and group and immunity test level for each emission and test standard specified by IEC60601-1-2 are summarized on EMC Declaration.
3. Accessories may affect EMC performance.

1.4.7 Safety Caution Statement

1. An identifying marking provided referring to instructions in IFU for battery pack intended to be changed only by service personnel using a tool.
2. A warning provided indicating replacement of Li-ion Battery pack IF incorrect replacement would result in an unacceptable risk.
3. A warning indicating that replacement by inadequately trained personnel could result in hazard.
4. Do not modify this equipment without authorization of the manufacturer.

Standards and Guidelines

The device meets the Safety and EMC requirements of the following standards and guidelines:

EMC and Safety Certified Standards for Powered Suction Pump

IEC 60601-1:2005 + A1:2012 + A2:2020

EN 60601-1:2006 + A1:2013 + A2:2021

IEC 60601-1-2:2014 + A1:2020

EN 60601-1-2:2015 + A1:2021

ANSI/AAMI ES 60601-1: 2005 & A1:2012 & A2:2021

CAN/CSA C22.2 No.60601-1: 14+A2:22(R2022)

Certified by NRTL_MET Classified

AS/NZS IEC 60601.1:2015 + Amd 1:2022

IEC 60601-1-11:2015 + A1:2020

EN 60601-1-11:2015 + A1:2021

EN ISO 14971:2019

Safety Certified Standard for Secondary Li-Ion Battery Pack

IEC/EN 62133-2:2017 and UN38.3_v7

1 Important Notes for Safe Use

Warranty terms and conditions

1.5 Warranty

The manufacturer guarantees the safety and correct functioning of the Powered Suction Pump only under the following conditions:

- The device is used for the intended purpose and maintained only in accordance with the information provided by these Instructions for Use.
- Only original spare parts or accessories approved by the Manufacturer are used.
- No structural changes are made to the device.
- Inspections and maintenance work are carried out by certified personnel according to specified time intervals.

1.6 Safety Information and Symbols of These Instructions for Use

The safety instructions in this Instructions for Use are marked with symbols and key words. Signal words like WARNING, CAUTION or ATTENTION designate the classification of the risk.

1.6.1 Identification of Risks of Injury

WARNING

Means a hazardous situation, which may lead to death or severe injury.

CAUTION

Means a hazardous situation, which may lead to minor or severe injury.

1.6.2 Identification of Material Damage

ATTENTION!

Describes a situation that could lead to property damage.

1.6.3 Identification of Additional Information

NOTE!

Means application tips and useful information.

Warning of damage to equipment surfaces

Advice to wear safety goggles

Advice to wear safety gloves


Advice to wear mouth and nose protection.

1 Important Notes for Safe Use

1.6.4 Additional Symbols to the Safety Information

Additional symbols to the safety information are those listed below

	Manufacturer		Country and date of manufacture
	Authorized Representative in the European Community		Unique Device Identifier
	Catalogue Number		US NRTL - MET Classified Mark
	Serial Number		Medical Device
	Follow Instructions for Use		Caution (ISO 7000-0434A)
	The CE mark indicates compliance with European harmonized legislation, and the numbers represent the supervision of a Notified Body		Type BF Applied Part
	The CE mark indicates compliance with European harmonized legislation		Waste Electrical and Electronic Equipment (WEEE Logo)
	Double Insulated, Class II AC Adapter		Temperature limit
IP22	Protected against ingress of solid foreign objects ≥12.5mm diameter. Protected from water spray less than 15 degrees from vertical		MR unsafe-Keep away the device from Magnetic resonance imaging (MRI) equipment
			Packaging unit
	Federal (US) law restricts this device to sale by or on the order of a physician.		Consult instruction for use
	Model Number		

1.7 Safety Instructions Before Use

- (1) The Powered Suction Pump can be administered only by persons or qualified medical staff who have been trained in its operation according to the instruction guidelines issued by the manufacturer.
- (2) Before using the Powered Suction Pump as a vacuum source and treatment system, please read the indications, warnings, precautions, and contraindications.
- (3) Check function of the unit prior to each use.
- (4) Never connect the power supply adapter to defective power sockets.
- (5) Keep power supply adapter and cable away from external heat sources. DO NOT cover the power supply adapter.

1 Important Notes for Safe Use

- (6) The device should not be charged or started up:
 - If the power cord or plug are defective;
 - If the device is not functioning properly;
 - If the device has been damaged/dropped;
 - If the device has been dropped into water;
 - If obvious defects might restrict safe operation
- (7) The Powered Suction Pump must be placed carefully and securely at the patient's bedside with optional VT accessories. An optional carrying bag is available for mobile use; however, it is the responsibility of the clinician or trained caregiver to determine if the patient's condition allows for mobile use.
- (8) It must be ensured that in between different patients use the carrying bag is disinfected or a new carrying bag is used.
- (9) The Powered Suction Pump must never be used to remove explosive gases and inflammable or corrosive fluids.
- (10) The unit must not be operated in damp rooms or when taking a bath or shower.
- (11) Avoid moisture on plug and switches. Never plunge the unit into water or liquids, not even when it is switched off.
- (12) The unit must not be operated in splash water range or in locations where there is a danger of explosion.
- (13) Operation of the Powered Suction Pump is possible while the battery is charging.
- (14) Pay attention to the ambient conditions described in the technical data. If the unit is operated at ambient temperatures outside the stated temperature range (see "Technical Data"), the performance may be reduced and the unit or the electronics and battery may get damaged.
- (15) The unit should be operated on a firm, level surface.
- (16) When the device is switched on, DO NOT leave it unattended.
- (17) Parts of the unit should be checked for correct function and safety-related defects at regular intervals. Please refer to the service manual.
- (18) The Powered Suction Pump must be switched off and disconnected from the power supply adapter before cleaning and maintaining unit.
- (19) The Powered Suction Pump is a medical device; it is not a toy. Keep away from children and pets as they can damage the dressing and therapy unit and affect performance.
- (20) Keep unit free of dust and lint.
- (21) Advise patient to NOT SUBMERGE therapy unit or dressing in liquid and to ensure therapy unit is not pulled into a tub or sink where it may become submerged.
- (22) Only use original, genuine Carilex Medical accessories and spare parts.
- (23) The unit must be used only with a genuine collection canister.

CAUTION

AVOID ELECTRIC SHOCK, DO NOT OPEN SUCTION PUMP!

WARNING

Electric shock!

The touching live parts can result in a burn by an electric shock. Check for damage of the plug and the main power cable of the power unit before connecting.

ATTENTION

Device fall may bump on the patients causing swelling and pain, please fix the device on the I.V. stand or bed frame, carrying bag.

2 Delivery and Storage

2.1 Packaging

The Powered Suction Pump is supplied with sturdy cardboard packaging. All packaging materials are recyclable and can be separated:
Packing: Cardboard, waste paper.

2.2 Models

Catalogue number	Model	Spec.
S1002-0022	VT · 200-i	Good for general NPWT purpose
S1002-3022	VT · 200-i	VT · 200-i P.U.

2.3 Delivery Control

Check immediately after delivery of the device :

- The completeness of the delivery
- The delivery status of the device

The Powered Suction Pump is delivered with the following components:

1. One Powered Suction Pump control unit
2. One 300ml canister
3. External AC Adapter
4. Instructions for Use


If the delivery is incomplete or the device and/or the packaging are damaged, in particular that caused by moisture or water, promptly inform the carrier as well as the distributor.

To fully charge the battery, and prior to the first start, attach the Powered Suction Pump to a wall outlet.

2.4 Operating & Transport & Storage Conditions

Recommended environmental conditions:

- For Operating Conditions
Temperature range: 5°C (41°F) to 35°C (95°F)
Relative Humidity range: 15% to 90%
Atmospheric Pressure range: 700 hPa to 1060 hPa
- For Transport and Storage Conditions
Temperature range: -25°C (-13°F) to 70°C (158°F)
Relative Humidity range: 0% to 90%

For long-term storage the power unit should be covered with a dust protector and the battery needs to be recharged every three months.

⚠ ATTENTION

Storage of the Powered Suction Pump

- Keep away from high voltage
- Keep away from humidity
- Keep away from heat
- Keep out of the reach of children
- Properly store in its box
- Do not store with other equipment


3 Device and Functional Description

3.1 Device Description

Powered unit

The power unit **1** is used as the housing for the compressed air unit as well as the compressed air system. It features:

- Display panel **2** for information.
- A control panel **3** with buttons to turn unit on/off and select functions.
- Canister **4** to collect fluid.

Power cord and charge **5**

Select the correct AC adapter plug for your country then plug into the main outlet.

ATTENTION

When changing to a different AC adapter plug, you will hear a “click” sound when the plug is firmly connected. Failure to fully connect will cause power failure.

NOTE


Please read the instruction for use prior to first use.


P/N R-NA-2(R)
North America
China
Japan

P/N R-SAA-2(R)
Australia

P/N R-UK-2(R)
United Kingdom
Hong Kong
Singapore

P/N R-EU-2(R)
Europe
South America


3 Device and Functional Description

3.2 Functional Description

Carilex Powered Suction Pump is a Negative Pressure Wound Therapy device that has been prescribed by your healthcare provider. This device has shown that it may help promote the healing of several different kind of wounds. When in use, negative pressure (suction) is delivered to the wound through the pump. The suction of the pump will help remove excess fluids from the wound. A special dressing will be placed onto your wound by healthcare professional and tube will be connected from your wound to the canister on the pump. After the dressing and tube are correctly applied and connected, turn on the Powered Suction Pump and set to the pressure setting that is prescribed by your healthcare provider. The canister will then collect excess fluid.

3.2.1 Panel of the Power Unit

- (a)  On / Off
- (b)  Arrow up
- (c)  OK / Unlock / Stop
- (d)  Mute / Return
- (e)  Display Screen


3.2.2 Function of the Power Unit

The vacuum air suction pump in the power unit sucks the air from the wound through the connecting tube and dressing to create a negative pressure environment in the wound.

3.2.3 Therapy Mode


Continuous mode:

When employed via the screen panel, the therapy unit will apply negative continuously. Default pressure for continuous mode is -125mmHg.


Intermittent mode:

When employed via the screen panel, the therapy unit changes between a high and low pressure range in a fixed time interval. Default pressure is -80mmHg 5 minutes / -20mmHg 2 minutes.


Instill mode:

This device allows for the optional use of irrigation solution to the wound bed. Irrigation solutions should only be used with direct supervision of a trained clinician. Default instill pressure is -20mmHg 10 minutes/ continuous -125mmHg 1 hour in one cycle.

4 Getting Started

3.2.4 Symbols on the Display

	Battery status
	Low power
	Battery cannot charge. Please contact your distributor for service.
	Charging
	Panel locked
	Error indicator
	The machine is running
	Call for service
	Battery malfunction but still working

4 Getting Started

NOTE

Before using, inspect the dressing kits to make sure there is no damage to the packaging, which may compromise the sterility of the contents. DO NOT use the contents of a damaged package; instead, dispose of it properly.

4.1 Battery

Charging instruction for battery

Select the correct AC adapter plug for your country. Connect the power supply adapter to the AC main outlet and open the protective rubber cover of the unit for Direct Current (DC) socket. Bars on the battery display will indicate that charging is occurring. Long battery life lasts up to 24 hours, for uninterrupted therapy and portability.

ATTENTION

When changing to a different AC adapter plug, make sure you hear a “click” sound to confirm that the plug is firmly connected. Failure to do so will cause power failure. Upon initial receipt of the Powered Suction Pump and prior to first use, charge the battery for at least 6~8 hours in order to create optimum battery memory and maximize the number of charging cycles over the lifetime of the battery.

4 Getting Started

ATTENTION

The battery must be fully charged prior to first use of the Powered Suction Pump. If therapy unit is in warehouse / inventory and not used for more than three months, the battery needs to be recharged. Powered Suction Pump is equipped with a Li-ion battery. The battery will discharge depending on the run time of the therapy unit and through extended periods of inactivity.

Depending on usage, the battery life cycle is stated as 7500 hours. Storage and usage of the battery must be within the temperature ranges stated in the section 2.4.

ATTENTION

The pump can still be charged when the battery is broken or unconnected for urgent needs, but users are not suggested to use the pump under these circumstances.

NOTE

For long term storage, disconnect the power adapter from the electrical outlet. When the machine is not used, it is recommended to fully charge the battery every three months to maintain battery life and ensure proper functioning.

NOTE

Dispose of the battery according to local or facility guidelines.

To protect the environment, dispose of useless devices at appropriate collection sites according to national or local regulations.

4.2 Collection Canister

Always make sure the canister is properly inserted. You will hear a "click" for proper installation, and it must remain in an upright (display side up) position during use.

The Powered Suction Pump is protected against penetration from solid / fluid substances by a hydrophobic membrane integrated with an activated carbon filter. If this filter fails, the Powered Suction Pump must be replaced.

The VT · 200-i system is designed for detecting when canister is full. When the liquid absorber reaches the canister full level, audible and visual indicators will be triggered and the message indicator "canister full" will be seen on the display panel. To avoid breaking of the suction inlet on the canister, do not pull the tubing of the canister horizontally.

The collection canister is to be properly discarded when full; it must be replaced after every patient use. The canister should also be inspected and replaced weekly or between patient users or otherwise as needed.

4 Getting Started

NOTE

The canister is designed for single patient use only. DO NOT re-use the canister to avoid the cross infection in between patients.

Safety Notice :

- Keep children clear from the extra-long cable and hoses to avoid entanglement and strangulation.
- Choking hazard caused by small parts being inhaled or swallowed. Identify any loose or detached small parts and keep them away from children.
- Identify any rubber or latex potential allergic reactions to materials used in the equipment.
- Contact Injuries-Check for any skin irritation due to prolonged contact with the equipment.
- Protection against strangulation or asphyxiation. A medical professional shall provide the means to control the risk of strangulation and asphyxiation of the patient and others by routing wires appropriately.

4.3 Dressing

Only the VT Dressing Kit is to be used in conjunction with the VT · 200-i system and must be in sterile condition. VT Dressing Kit should be applied in accordance with the Dressing Kit Instructions for Use, supplied with the dressings.

A physician or trained caretaker should perform an intensive, thorough wound cleaning prior to applying the dressing.

Perform routine dressing checks and changes every 48 hours or according to the facility protocol or physician's order.

Before using, inspect the dressing kits to make sure there is no damage to the packaging, which may comprise the sterility of the contents.

CAUTION

DO NOT use the contents of a damaged package; instead, dispose of it properly.

DO NOT re-use the dressing kit to avoid cross infection. Follow IFU for dressing kits.

Wound infection that may increase the risk associated with wound worsen, sepsis or osteomyelitis. Please be careful to seal the dressing, or, contact your caregiver when the dressing come off from the skin.

4 Getting Started

4.4 Lay Operator Briefing Information

The equipment and accompanying Operation Manual should be simple to understand and straightforward to use for the Lay Operator, according to IEC60601-1-11 regulation, a Lay Operator needs to have at least 8 years of education. The healthcare professional should brief the Lay Operator on the use of the equipment and any precautions to be taken, including:

- Precautions to be taken in the event of changes in the performance of the equipment.
- Precautions to be taken regarding the exposure of the equipment to reasonably foreseeable environmental conditions.(e.g. magnetic and electromagnetic fields, external electrical influences, electrostatic discharge, variations in pressure etc.)
- Information about medicinal substances that the equipment is designed for, including any that the equipment is NOT designed for.
- Information about medicinal substances or blood products incorporated into the equipment as an essential part.
- The accuracy of equipment with a measuring function.

carilex®

5 Operation Procedure

"Carilex" Canister lists		
1	S1001-6020	Canister for VT· 100 - 500ml lock
2	S1001-6040	Canister for VT· 100 - 1000ml lock
3	S1001-6060	Canister for VT· 100 - 300ml lock
4	S1001-6100	Canister for VT· 100 - 300ml luer lock
5	S1001-6110	Canister for VT· 100 - 500ml luer lock

"Carilex" Dressing lists		
1	S1001-2040	VT Dressing Kit S Carilex
2	S1001-2050	VT Dressing Kit M Carilex
3	S1001-2060	VT Dressing Kit L Carilex
4	S1001-2100	VT Spiral Dressing Kit Black Lock
5	S1001-2120	VT Dressing Kit XL Carilex
6	S1001-2130	VT Drape and Port Set Kit S Lock
7	S1001-2140	VT Drape and Port Set Kit M Lock
8	S1001-2150	VT Drape and Port Set Kit L Lock
9	S1001-2200	VT Dressing Kit TS S Carilex
10	S1001-2210	VT Dressing Kit TS M Carilex
11	S1001-2220	VT Dressing Kit TS L Carilex

5.1 Check Points before Using Powered Suction Pump

Before using the Powered Suction Pump, it is important to check for the following:

- Damage to the power cord and plug
- Damage to the pump
- Complete of the packaging
- Battery status

5.2 Insert Canister into the Powered Suction Pump

(Please follow the instructions on 4.2 Collection Canister section)

Insert the canister into the Powered Suction Pump connecting port.

Make sure the level marker of the canister on the Powered Suction Pump is on the same side as the display screen of the pump.

5.3 Connect the Powered Suction Pump Tubing Connector to the Carilex VT Dressing Kit Tubing Connector, and Screw Tightly to Ensure a Tight Airlock.

WARNING

Electric Shock!

Touching live parts can result in a death or serious injury by an electric shock. Check for damage to the plug and the main power cable before connecting.

- Do not use damaged components.

5 Operation Procedure

5.4 Using the Function Keys

1. Power 

Press  and hold button for 3 seconds to turn on the power unit. The LED will display a "green" color.
2. Mute / Return 

After receiving an audible notification, you may press  button it turn off. Press the button system to return to the previous menu. In the sub-menu, press  button to return to the main menu.
3. Arrow up 

In the continuous and intermittent setting mode, press  to select the desired pressure or minutes. -5mmHg or 1minute per step.
4. OK / Unlock / Stop 

Press  button to select item. In main-menu, press  button to get to the sub-menu of options. In therapy mode, press  to stop the therapy.

When the control panel is locked, press  button for 3 seconds to unlock.


5.5 Turn the Unit On / Off

1. Turn on the unit


To turn on the unit, press the power button on the control panel for 3 seconds until the screen is displayed as Figure (a).
2. Turn off the unit

To turn off the unit, press the power button on the control panel for 3 seconds. If the therapy unit is locked, press  for 3 seconds to unlock it, then press the power button for 3 seconds to turn off the unit. The display screen will go blank after the unit is turned off.


(a)

(b)

(c)

5.6 Therapy Mode Setting

1. To select therapy mode, switch on the therapy unit with the power key.
2. When the screen shows VT · 200-i, immediately press  key if you want to continue with the previous setting, or wait for 5 seconds until the screen displays Figure (b).
3. Use  keys to move the cursor to the prescribed therapy mode and press  to confirm.

Continuous mode setting (c)


1. Press  to change the constant pressure level from -30mmHg to -200mmHg. The default pressure for continuous mode is -125mmHg.
2. Press  to start the therapy.

5 Operation Procedure


(d)

(e)


(f)


(g)

(h)

(i)

Intermittent mode setting (d-e)

1. Press to change the high pressure level from -50mmHg to -200mmHg. Default pressure for high is -80mmHg. Press to confirm the high pressure setting.
2. Press to change the high pressure therapy time from 1 to 10 minutes. Default time for high pressure is 5 minutes. Press to confirm the high pressure time setting.
3. Press to change the low pressure level from -20mmHg to -80mmHg. Default pressure for low pressure is -20mmHg. Press to confirm the low pressure setting.
4. Press to change the low pressure therapy time from 1 to 10 minutes. Default time for low pressure is 2 minutes. Press to confirm the low pressure time setting.

Instill mode setting (f-i)

1. Press to select continuous pressure level from -80mmHg to -200mmHg. Default pressure for therapy is -125mmHg. Press to confirm the therapy pressure.
2. Press to select therapy time (1,2,4,6,8,10,12 hours). Default time for therapy is 1 hour. Press to confirm the therapy time.
3. Press to select instill therapy pressure from -20mmHg to -80mmHg. Default pressure for instill is -20mmHg. Press to confirm the instill pressure.
4. Press to select instill time (10,20,30,40,50,60 minutes). Default time for instill is 10 minutes. Press to confirm the instill time.
5. Press to select instill cycle time from 1 to 10 times. Default cycle time is 1. Press to confirm the instill cycle time.
6. Press to select which therapy cycle to start with.

5 Operation Procedure

5.7 Starting Negative Pressure Wound Therapy (NPWT)

The "Carilex" VT · 200-i model allows for the optional use of wound irrigation solution to the wound bed, cleanse with the wound irrigation solution of topical wound. The Carilex VT Dressing Foams are hydrophobic properties, provide improved fluid distribution within, and removal of fluids, including wound exudate, irrigation fluids, body fluids, infectious material and tissue debris which may promote wound healing.

**NOTE**

Wound irrigation solution should only be used with direct supervision of a trained clinician. Wound irrigation solution is not intended for use in a homecare setting.

1. Prior to dressing application, thoroughly clean the wound and dry with aseptic techniques, debride any necrotic or devitalized tissue in accordance with your institutions protocols.
2. Cut the supplied foam to fit the wound. Place the foam into the wound.
3. Apply the drape over the foam allowing at least 5 cm margin to the wound edge. Push down around the drape to seal against the patient's skin.
4. Pull a small amount of drape up from the center of the wound location and cut with scissors to leave a hole in the drape of approximately 2 cm.
5. Apply to the Port and seal to the drape.
6. Connect the VT Dressing tubing to the VT therapy unit Canister tubing.


Continuous or Intermittent mode, please skip to step 12.

Select the Instill mode, please see the following steps

7. Obtain physician order for instill therapy mode and delivery drip rate.
8. Make sure the instill fluid supply remains off until therapy is started and target pressure achieved.
9. Select desired location for instill port. Cut another hole of approximately 2cm with scissors in the drape, which the instill port would be placed.


**NOTE**

Please mind the two ports need to keep a certain distance.

**NOTE**

The drip rate does not need to be exact with continuous wound instill.

5 Operation Procedure


10. Connect the VT instill tubing to the IV port (Male luer lock \varnothing 3.9mm) tubing. Close irrigation solution delivery set clip completely.
11. Hang irrigation solution bag on IV pole higher than the wound.
12. Check tubing connectors to ensure they are fully engaged and locked.
13. Ensure the tube clip are released from the tubing.
14. Follow these Instructions for Use of the Carilex VT Instill Therapy unit before activate the device.
15. Turn device ON and allow dressing to reach set pressure.
16. Check the seal of the dressing once therapy is activated.

NOTE

Larger dressings will take longer to evacuate.



17. Open the clip and allow irrigation solution to flow until the irrigation solution begins to move through the suction tubing and into the canister.

18. Set the drip rate per physician order.
19. Lock the IV tube clip after the instill mode is finished with audible beeps. If instill cycle is selected more than 1 time, please repeat the above steps.

5 Operation Procedure

(j)

Others

1. When left untouched for 3 minutes, the panel will lock. To unlock the panel, press  for 3 seconds.
2. If the therapy needs to be stopped, press  to stop. Press  again to restart.
3. If the prescribed pressure / therapy mode needs to be changed, press  to unlock the panel, and then press  to go back to the original setting.
4. During the therapy, press  and the setting details (j) will be displayed for 5 seconds.
5. To save power, the screen will turn dark after 3 minutes.

ATTENTION

If the suction that regulates the pressure does not perform properly, the device may apply higher pressure than the intended value and may cause bleeding and blood loss. If the device stops working or pressure is too low, it can cause delayed healing, loss of viable tissue or contamination of the wound.

ATTENTION

If the suction that regulates the pressure does not perform properly, the device may apply higher pressure than the intended value and may cause loss of soak functionality.

ATTENTION

If the suction that regulates the pressure does not perform properly, the device may apply lower pressure than the intended value and may cause over soak to induce contamination of the wound.

5 Operation Procedure

5.8 Mute / Return

Press  to mute the pump when acoustic signals occur.

When the pump is operating in CONTINUOUS MODE, press  to adjust the pressure level.

When the pump is operating in INTERMITTENT MODE, press  to adjust the pressure level.

Press  two times to go back to therapy mode selection menu.

The Powered Suction Pump is made from various electronics and plastics. When the Powered Suction Pump is ready for disposal, follow local governing guidelines regarding appropriate and proper disposal procedures of the device components.

The used canisters, tubes and dressings should be disposed according to the local or facility guidelines for handling infected or bio- hazardous materials. None of the items should be disposed or co-mingled with household or facility refuse. Incorrect disposal may have harmful effects to the environment and public health.

5.9 Change Powered Suction Pump Canister

The canister has to be changed by visual check or according to the instructions on the display.

The VT · 200-i system is designed for detecting when the canister is full. When the liquid absorber reaches the canister full level, sound and visual indicators will trigger and the message "canister full" will show on the display panel.

Do not pull the canister tubing horizontally to avoid the breaking the suction inlet on the canister.

Properly discard the collection canister when full; it must be replaced after every patient use. Also, inspect and replace the canister weekly and/or between patient uses, or otherwise as needed.

ATTENTION

Avoid cross infection between patients and DO NOT re-use the canister.


The procedure of changing canister:

1. Wash hands and wear disposable gloves.
2. Switch the Powered Suction Pump off by pressing  for 3 seconds.
3. Close the clip on the tube, both the dressing and canister's side, and then disconnect the canister connector from the dressing.
4. Detach the full canister. Discard the used canister by following local governing guidelines.
5. Attach a new canister to Powered Suction Pump. (Refer to 4.2 Collection Canister section for instructions.)

6 Error Indicators

6.1 Error Indicators

If the Powered Suction Pump detects any of the following Error Messages, the display screen turns on with an acoustic warning signal simultaneously. Press the mute key to turn off the acoustic indicator. The orange light will blink (1 sec on/ 1 sec off) until the issue is resolved.

Error Message	Display	Possible causes	Remedy
Leakage		Dressing is not tight. Tube is not well connected or leakage occurred in the dressing.	Press  for 3 seconds to turn off the power suction pump. Check the system for leakage. Turn on the power suction pump again after leakage issue is resolved.
Blockage		Tubing is kinked or clip closed. Tube clogged.	Press  to continue therapy. Check the system for blockages. If the blockage issue is still not resolved, blockage indicator will happen again. 
Canister Full		Canister full	Press  for 3 seconds to turn off the power suction pump. Replace with new canister.
Low Battery		Battery low	Press  for 3 seconds to mute the acoustic signal. The remaining battery time is approximately 10 minutes. Charge the battery. The visual symbol blinks and the acoustic signal will be activated until the battery is empty.
Call for Service		Failure of PCBA, Battery, or Vacuum Motor	Press  for 3 seconds to turn off the power unit. Contact your authorized local distributor for assistance.

7 Application of the Negative Pressure Wound Therapy

7.1 Application of the Negative Pressure Wound Therapy

Comply with all hygiene regulations!

Basic cleaning

The components of the Negative Pressure Wound Therapy are not supplied in a sterile condition. Clean and disinfect the components before the first use and between each patient.

If the Powered Suction Pump fails to function normally for one hour, regardless of the circumstances, stop using the Powered Suction Pump immediately. Turn off the power, remove the canister, and remove the VT dressing from the wound site and replace it with traditional wound dressing.

Do not leave the idle Powered Suction Pump with dressing in the wound site for more than one hour to avoid risk of cross contamination.

Check the dressing checked regularly to avoid blockage and leakage.

A loose power cord may cause tripping and serious injury. Make sure any cords and tubing is safely stowed.

Improper use of the Powered Suction Pump may cause pain and injury to the patient. Excessive negative pressure or an infection of the wound may cause pain and injury to the patient.

8 Disinfection and Cleaning

In order to prevent cross-contamination, the disinfection and cleaning of the entire Negative Pressure Wound Therapy unit must be completed before first use and between patients. If there is a identifiable disease according to the Federal Law concerning epidemics, consult a hygiene expert prior to disinfection and cleaning.

Standard Precautions are designed to reduce the risk of cross contamination of microorganisms from both known and unknown sources. Regardless of the patients' diagnosis or presumed infection status, this Precautions should be suitable to apply to all of them. Especially when users who had contact with blood and body fluids. This also includes secretions and excretions (except sweat) regardless of whether blood is visible or not, non-intact skin (i.e., open wounds) and mucous membranes.

All the disposable items such as tubing, connectors, clamps, used canisters, used dressing..etc. should discard in accordance with local medical waste disposal regulations. Improper disposal may run the risk of regulatory non-compliance.

WARNING

Electric shock!

Water has a high electrical conductivity. Contact with liquid under voltage can lead to a fatal electric shock. For the disinfection and cleaning operations:

- Turn off the power unit.
- Unplug it from the power socket.

CAUTION

Health hazard!

The contact with contaminated cleaning fluids can cause infections. Disinfectants can contain harmful substances.

Please follow the Instructions for Use of the manufacturer of the disinfectant and the hygiene of the operator during the disinfection and cleaning. Wear personal protective equipment:

- Safety glasses.
- Protective gloves.
- Mouth and nose protective.

ATTENTION

Incompatible cleaning agents!

Parts of the Negative Pressure Wound Therapy device are made of plastic. Solvents can damage plastic and coatings. Strong acids or alkalis can cause them to become brittle.

8 Disinfection and Cleaning

Hygiene requirements of the operator

8.1 Disinfection and Cleaning

Notify the operator about the measures which apply to Negative Pressure Wound Therapy and the hygiene directives for disinfection. The disinfection of the Negative Pressure Wound Therapy device or parts of it can be performed only by cleaning experts who are familiar with the hygiene requirements of the institution.

Disinfection procedures

Please follow the procedure required by your local health authority.

Disinfection procedures

Disinfection operation

Manual disinfection by wiping is carried out in three steps:

- Pre-disinfection
- Cleaning
- Controls

Pre-disinfection

- (1) Wear surgical gloves and surgical mask.
- (2) Wash hands with 75% alcohol or other cleansers in accordance with local Competent Health Authority regulation.
- (3) Wipe the surfaces with disinfection.
- (4) Allow the disinfectant to take effect according to the manufacturer's instructions.

Cleaning

- (1) Use 75% alcohol to wipe off dirt and dust accumulations for disinfection.
- (2) Wipe the surfaces with a clean soft cloth and clear water.
- (3) Dry all the surfaces with a clean soft cloth.

Controls

Check the function of the power unit.

CAUTION

For repair, contact your local distributor.
Please follow hygiene control regulations of your local government authority.

9 Care and Maintenance

9.1 Inspection

The operator must check the condition of the Negative Pressure Wound Therapy for each use; during use by patients; and at least once in a year especially in relation to the following:

- Function of the power unit with pressure control and canister control markings.
- Condition of the compressed air hoses and connections.

9.2 Maintenance

When the Powered Suction Pump is not in use, please recharge it every three months.

carilex®

10 Troubleshooting

Problem	Inspection Procedure	Possible Solution
1.) Power unit does not function	1.1) Check if power cord is firmly plugged into wall outlet	1.1) Secure power cord into wall socket
	1.2) Check if battery is empty	1.2) Connect the power supply adapter to the electrical outlet to recharge battery
	1.3) Check if the unit is in the Setting mode.	1.3) Complete the settings
2.) Insufficient performance	2.1) Check if the tube is partially kinked or the clip is engaged	2.1) Release clip, or remove kinks or check the system for blockages.
	2.2) Check system for leaks	2.2) Connect the tubing / canister properly. Seal dressings properly.
	2.3) Check if battery is almost empty	2.3) Charge the battery
	2.4) Check that the height of the suction unit is not more than one meter above the wound	2.4) Move unit to within one meter height of wound
3.) No suction	3.1) Check if tubing is blocked	3.1) Release clip, remove kinks or check the system for blockages. Remove clogs by applying a new dressing kit or pressing  &  to start the suction at full speed to clear the blockage. ⚠ The blockage removal function is designed to aid the removal of blockage. Due to the varying condition of individual wounds, this function may not guarantee complete and successful removal of any blockage. ⚠ The blockage removal function will create high pressure suction, which may cause excessive bleeding. Please consult with the physician prior to activating the blockage removal function.
		3.2) Check if the canister is full

If the troubleshooting procedures do not return the system to normal performance, stop using the system immediately and contact the authorized distributor for technical service.

11 Technical Data

	<p>The VT · 200-i is suitable for continuous operation.</p> <p>ModelVT · 200-i</p> <p>Dimensions (WxHxD).....18 x 17.9 x 8.9cm</p> <p>Weight (with empty canister)1.35kg</p> <p>Negative Pressure Range..... Continuous: -30mmHg~-200mmHg Intermittent: High -50mmHg~-200mmHg Low -20mmHg~-80mmHg Instill: -20mmHg~-80mmHg</p> <p>Suction Flow Rate.....4~5.5L/min</p> <p>Power Consumption.....30W max.</p> <p>Operating Noise.....50dBA at 190cm</p> <p>Nominal voltage9.1Vdc</p> <p>Medical device....Class IIb (Comply with 2017/745 as amended by 2019/5) Class II (Comply with 21 CFR 878.4780)</p> <p>IP CodeIP22</p> <p>Applied Part.....Type BF</p>
Powered Suction pump	
	<p>ANSI/AAMI ES60601-1_v3.2</p> <p>IEC/EN 60601-1_v3.2</p> <p>Certified ModelEljintek_GMPU30UI-2 or GlobTek_GTM91120-3010.5-1.4-T2</p> <p>Electrical ClassClass II</p> <p>Power input.....100-240Vac 50/60Hz 1.0-0.6A</p> <p>Power output.....9.1Vdc 3.29A.</p>
AC Adapter	
	<p>IEC/EN 62133-2 : 2017</p> <p>Certified Model59H-0186-506-G or 59H-0186-508-G;</p> <p>Rechargeable Li-Ion Battery Pack; 2S2P</p> <p>Normal Rating.....7.2Vdc 5.4Ah or 7.2Vdc 6.6Ah</p> <p>Battery Using Time.....24 hours max.</p> <p>Charging Time.....8 hours</p> <p>Be aware of country-specific, local and facility regulations regarding disposal, especially of used batteries, electronics, and medical waste.</p> <p>Operating Conditions</p> <p>Temperature Range: 5°C to 35°C</p> <p>Relative Humidity Range: 15% to 90%</p> <p>Atmosphere Range: 700hPa to 1060hPa</p> <p>This system is not AP / APG protected.</p> <p>The NPWT must be decontaminated before disposal.</p> <p>Disposal of old electrical and electronic equipment - valid in the European Union: WEEE Directive 2012/19/EU.</p>
Li-Ion BP	

This symbol on the product or on its packaging indicates that this product should not be treated as household waste. Instead, this product should be taken to the appropriate place of disposal for the recycling of electrical waste and electronic equipment.

12 EMC Declaration

Declaration of Conformity

For EN 60601-1-2 (4th Ed.)

Company Name:	Carilex Medical, Inc.
Company Address:	No. 77, Keji 1st Road, Guishan District, Taoyuan City 33383, Taiwan (R.O.C.)
Product Name:	NPWT, Powered Suction Pump
Model No.:	VT · 100, VT · 200, VT · 200-i
Report Number:	ETC 22-08-RBO-034
Power Supply:	1) 100-240Vac 50-60Hz or 50/60Hz...External AC Adapter 2) 7.2Vdc.....Internal Li-Ion Battery Pack

Recommended separation distances between portable and mobile RF communications equipment and the ME equipment			
The Powered Suction Pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Powered Suction Pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Powered Suction Pump as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.1	0.1	0.2
0.1	0.4	0.4	0.7
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	11.7	11.7	23.3

Declaration – electromagnetic emissions		
The Powered Suction Pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Powered Suction Pump should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Powered Suction Pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The Powered Suction Pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

12 EMC Declaration

Declaration – electromagnetic emissions and immunity – for EQUIPMENT and SYSTEMS that are use in the professional healthcare facility environment or in the home healthcare environment

The Powered Suction Pump declaration – electromagnetic immunity

The Powered Suction Pump system is intended for use in the electromagnetic environment specified below.
The customer or the user of the Powered Suction Pump system should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance	
Conducted RF IEC 61000-4-6	3 Vrms ; 6 Vrms 150 kHz to 80 MHz	3 Vrms ; 6 Vrms 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the EQUIPMENT or SYSTEM including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.	
Radiated RF IEC 61000-4-3	3 V/m ; 10V/m 80 MHz – 2.7 GHz 80%	10V/m 80 MHz – 2.7 GHz 80%		
Proximity fields from RF wireless Communications equipment IEC 61000-4-3	27 V/m	385 MHz	Interference may occur in the vicinity of equipment marked with the following symbol. 	
	28 V/m	450 MHz		
	9 V/m	710 MHz		28 V/m
		745 MHz		9 V/m
		780 MHz		710 MHz
	28 V/m	810 MHz		28 V/m
		870 MHz		810 MHz
		930 MHz		870 MHz
	28 V/m	1720 MHz		28 V/m
		1845 MHz		1720 MHz
1970 MHz		1845 MHz		
28 V/m	2450 MHz	28 V/m		
	5240 MHz	2450 MHz		
	5500 MHz	5240 MHz		
9 V/m	5785 MHz	9 V/m	5500 MHz	
			5785 MHz	

Declaration – electromagnetic immunity

The Powered Suction Pump system is intended for use in the electromagnetic environment specified below.
The customer or the user of the Powered Suction Pump system should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV ±1 kV differential mode ±2 kV common mode	±0.5 kV ±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U_T ; 0 , 5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycle Single phase: at 0°	0 % U_T ; 0 , 5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycle Single phase: at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EQUIPMENT or SYSTEM requires continued operation during power mains interruptions, it is recommended that the EQUIPMENT or SYSTEM be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

13 Recommendations for Li-Ion Battery Pack

The following represents a typical, but non-exhaustive, list of suggestions provided by the equipment manufacturer to end-users.

1. Do not short-circuit a cell or battery. Do not store cells or batteries haphazardly in a box or drawer where they may short-circuit each other or be short-circuited by conductive materials. Do not dismantle, open, or shred secondary cells or batteries.
2. Do not remove a cell or battery from its original packaging until required for use.
3. Do not subject cells or batteries to mechanical shock.
4. In the event of cell leakage, do not allow the liquid to come into contact with skin or eyes. If contact is made, wash the affected area with copious amounts of water and seek medical attention.
5. Do not mix cells of different manufacture, capacity, size or type within this medical device.
6. Seek medical advice immediately if a cell or battery has been swallowed.
7. Keep cells and batteries clean and dry.
8. Wipe the cell or battery terminals with a clean dry cloth if they become dirty.
9. Secondary cells and batteries need to be charged before use.
10. After extended periods of storage, it may be necessary to charge and discharge the cells or batteries several times to obtain maximum performance.
11. The Li-Ion Battery Pack should be charged at temperature between 10°C and 45°C
12. The Li-Ion Battery Pack should be discharged at temperature between 10°C and 60°C
13. Do not use any charger other than that specifically provided for use with the equipment.
14. Do not use any cell or battery which is not designed for use with the equipment.
15. Battery usage by children should be supervised.
16. Always purchase the battery recommended by the device manufacturer for the equipment.
17. Do not leave a battery on prolonged charge when not in use.
18. Use only the cell or battery in the application for which it was intended.
19. When possible, remove the battery from the equipment when not in use.
20. Dispose of properly.

WARNING

The risk of battery thermal runaway may causing patient burn. Therefore, user shall follow the info of IFU, proper maintain battery.

carilex[®]

carilex®

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"Carilex" VT-200-i

Negative Pressure Wound Therapy System

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Declaration – electromagnetic emissions		
The Powered Suction Pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Powered Suction Pump should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Powered Suction Pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	The Powered Suction Pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

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The Powered Suction Pump declaration – electromagnetic immunity					
The Powered Suction Pump system is intended for use in the electromagnetic environment specified below. The customer or the user of the Powered Suction Pump system should assure that it is used in such an environment.					
Immunity test	IEC 60601 test level		Compliance level		Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms ; 6 Vrms 150 kHz to 80 MHz		3 Vrms ; 6 Vrms 150 kHz to 80 MHz		Portable and mobile RF communications equipment should be used no closer to any part of the EQUIPMENT or SYSTEM including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m ; 10V/m 80 MHz – 2.7 GHz 80%		10V/m 80 MHz – 2.7 GHz 80%		
Proximity fields from RF wireless Communications equipment IEC 61000-4-3	27 V/m	385 MHz	27 V/m	385 MHz	Interference may occur in the vicinity of equipment marked with the following symbol. 
	28 V/m	450 MHz	28 V/m	450 MHz	
	9 V/m	710 MHz	9 V/m	710 MHz	
		745 MHz		745 MHz	
	28 V/m	780 MHz	28 V/m	780 MHz	
		810 MHz		810 MHz	
		870 MHz		870 MHz	
	28 V/m	930 MHz	28 V/m	930 MHz	
		1720 MHz		1720 MHz	
		1845 MHz		1845 MHz	
28 V/m	1970 MHz	28 V/m	1970 MHz		
	2450 MHz		2450 MHz		
9 V/m	5240 MHz	9 V/m	5240 MHz		
	5500 MHz		5500 MHz		
	5785 MHz		5785 MHz		

Declaration – electromagnetic immunity			
The Powered Suction Pump system is intended for use in the electromagnetic environment specified below. The customer or the user of the Powered Suction Pump system should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV ±1 kV differential mode ±2 kV common mode	±0.5 kV ±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U_T ; 0 , 5 cycle At 0° , 45° , 90° , 135° , 180° , 225° , 270° and 315° 0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycle Single phase: at 0°	0 % U_T ; 0 , 5 cycle At 0° , 45° , 90° , 135° , 180° , 225° , 270° and 315° 0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycle Single phase: at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EQUIPMENT or SYSTEM requires continued operation during power mains interruptions, it is recommended that the EQUIPMENT or SYSTEM be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

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VORSICHT

For repair, contact your local distributor. Please follow hygiene control regulations of your local government authority.

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12 EMV-Erklärung

**Declaration – electromagnetic emissions and immunity –
for EQUIPMENT and SYSTEMS that are use in the professional healthcare facility
environment or in the home healthcare environment**

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Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance	
Conducted RF IEC 61000-4-6	3 Vrms ; 6 Vrms 150 kHz to 80 MHz	3 Vrms ; 6 Vrms 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the EQUIPMENT or SYSTEM including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.	
Radiated RF IEC 61000-4-3	3 V/m ; 10V/m 80 MHz – 2.7 GHz 80%	10V/m 80 MHz – 2.7 GHz 80%		
Proximity fields from RF wireless Communications equipment IEC 61000-4-3	27 V/m	385 MHz	Interference may occur in the vicinity of equipment marked with the following symbol. 	
	28 V/m	450 MHz		
	9 V/m	710 MHz		28 V/m
		745 MHz		9 V/m
		780 MHz		710 MHz
	28 V/m	810 MHz		28 V/m
		870 MHz		810 MHz
		930 MHz		870 MHz
	28 V/m	1720 MHz		28 V/m
		1845 MHz		1720 MHz
1970 MHz		1845 MHz		
28 V/m	2450 MHz	28 V/m		
	5240 MHz	2450 MHz		
	9 V/m	5240 MHz		
	5500 MHz	9 V/m		
	5785 MHz	5240 MHz		

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